




YS6010 Benchtop Grating Spectrophotometer

YS6010 spectrophotometer is a bench-top measurement instrument designed to evaluate the color and appearance of many measuring applications for objects' reflectance and also the transmittance of translucent or transparent liquids and solids. Equipping with 7 inches TFT capacitive touch screen display, full illuminants, reflective d/8 and transmissive d/0 geometr(including or excluding UV) makes YS6010 ideal for color analysis within R&D and laboratory environments.



7 inches TFT Capacitive Touch-screen



Store approximately 20000 test data



Wavelength range 360nm – 780nm



Built-in camera locating




PRODUCT FEATURES

1. Double Array 256 Image Element CMOS Sensor; Long life-span stable LED, UV LED and xenon lamp.
2. With reflectance and transmittance spectrum, accurate L^*a^*b value, good to calculate color formula and do precise color transmission.
3. Auto-identify measuring aperture. Freely switchable between 3 reflective measuring apertures: $\phi 25.4\text{mm}/8\text{mm}/4\text{mm}$. Users can also customize apertures.
4. Built-in temperature sensor to monitor and compensate the measuring temperature to ensure the measurement more precision.
5. Wavelength range 360nm – 780nm. Built-in 400nm cut off/420nm cut off/460nm cut off (only xenon lamp edition), more professional in UV measurement.
6. Independent light source detector, continuously monitor the change of light sources to ensure the light source reliable.
7. Multiple accessories, sample holders, fixation frame, suitable to more working condition.
8. More powerful extended functions with the PC software.



PRODUCT HIGHLIGHTS



Plastic Sample

Reflectance Measurement:
Adopt $D/8^\circ$ Geometry,
conform to ISO7724/1, CIE,
ASTM, DIN and JIS standard.



Transmissive Sample

Transmittance Measurement: Adopt
 $D/0^\circ$ Geometry, conform to ISO, CIE,
ASTM and DIN standard.



Liquid, Powder, Solid Sample Measurement

360-780nm wide wavelength
range measurement,
professional UV sample
measurement.



Measuring Apertures

Auto-identify the apertures,
 $25.4/8/4\text{mm}$ three apertures
freely switch. Special aperture
can be customised.

APPLICATION INDUSTRY

YS6010 benchtop spectrophotometer is designed for accurate analysis and transmission of laboratory color and appearance of solids, liquids, pastes, powders, pills, and granules in a more comprehensive, streamlined process. YS6010 is widely used in paints, inks, textiles, garments, printing and dyeing, printing etc industries for color transfer and quality control.



Automobile



Leather



Plastics




Paint



Food stuff



Laboratory



Others

YS6010 Specification

Optical Geometry: Reflective : d/8°(diffused illumination, 8-degree viewing angle);SCI&SCE / SPIN&SPEX;Including / excluding UV measurements
Transmissive : d/0°(diffused illumination, 0-degree viewing angle);SCI&SCE / SPIN&SPEX;Including / excluding UV measurements
Comply to standard CIE No.15, GB/T 3978, GB 2893, GB/T 18833, ISO7724/1, ASTM E1164, DIN5033 Teil7

Integrating Sphere Size:Φ154mm

Light Source Device:360~780nm Combined LED Lamp, 400nm Cut-off light source

Application:It is used for accurate analysis and transmission of laboratory color
Apply in paints, inks, textiles, garments, printing and dyeing, printing etc industries for color transfer and quality control

Spectral Separation Device:Concave grating

Detector:256 Image Element Double Array CMOS Image Sensor

Wavelength Range:360~780nm

Reflectance Range:0~200%

Wavelength Pitch:10nm

Measuring Aperture:Reflective : Φ30mm/Φ25.4mm, Φ10mm/Φ8mm, Φ6mm/Φ4mm;Transmissive : Φ30mm/Φ25mm;

Specular Component:Reflectance:SCI&SCE;Transmittance: SCI&SCE

Color Space:CIE LAB,XYZ,Yxy,LCh,CIE LUV,Musell,s-RGB,HunterLab,βxy,DIN Lab99

Other Colorimetric Index:WI (ASTM E313, CIE/ISO, AATCC, Hunter),YI (ASTM D1925, ASTM 313),TI (ASTM E313, CIE/ISO),MI (Metamerism Index),
Staining Fastness, Color Fastness, Color Strength,Opacity,8° Glossiness, Gardner Index, Pt-Co Index, 555 Index

Color Difference Formula: $\Delta E^*ab, \Delta E^*uv, \Delta E^*94, \Delta E^*cmc(2:1), \Delta E^*cmc(1:1), \Delta E^*00, \text{DIN}\Delta E99, \Delta E(\text{Hunter})$

Illuminants:D65,A,C,D50,D55,D75,F1,F2,F3,F4,F5,F6,F7,F8,F9,F10,F11,F12

Displayed Data:Spectrogram/Values, Chromaticity Values, Color Difference Values/Graph, Pass/Fail Result, Color Offset

Measurement time: About 3s (Measure both SCI & SCE about 6s)

Repeatability:Spectral reflectance: Φ25.4mm/SCI, Standard deviation within 0.05% (400 nm to 700 nm: within 0.04%)

Chromaticity value:Φ25.4mm/SCI, Standard deviation within ΔE^*ab 0.02 (When a white calibration plate is measured 30 times at 5 second intervals after white calibration)
Spectral Transmittance: Φ25.4mm/SCI, Standard deviation within 0.05% (400 nm to 700 nm: within 0.04%)

Chromaticity value:Φ25.4mm/SCI, Standard deviation within ΔE^*ab 0.03 (When a white calibration plate is measured 30 times at 5 second intervals after white calibration)

Specular Component: Reflectance:SCI&SCE;Transmittance: SCI&SCE

Observer Angle: 2°/10°

Inter-instrument Error:Φ25.4mm/SCI, Within ΔE^*ab 0.15(Average for 12 BCRA Series II color tiles)

Size:370×300×200mm(L*W*H)

Weight:9.6kg

Power :AC 24V, 2A Power adapter power supply

Illuminant Life Span:5 years, more than 3 million times measurements

Display:7-inch TFT color LCD, Capacitive Touch Screen

Data Port:USB, Print serial port

Data Storage Capacity:Standard: 2000 Pcs; Sample: 20000 Pcs(One PCS can include both SCI and SCE)

Language:English & Chinese

Working Environment:Temperature: 0~40°C; Humidity: 0~85% (No Condensation)

Storage Environment:Temperature: -20~50°C; Humidity: 0~85% (No Condensation)

Standard Accessory:White and Black Calibration Board, Checking Green Board, Sample Holder, Φ4mm, Φ8mm, Φ25.4mm Aperture,
Power Adapter, USB Cable, User Guide, PC Software

Optional Accessory:Micro-printer, Transmissive Test Component

Notes: The specifications are subject to change without notice



SHENZHEN ThreeNH TECHNOLOGY CO., LTD.

Address:F/6, Block 5B, Skyworth Inno Valley, Tangtuo 1st Road,
Shiyan, Baoan District, Shenzhen, P.R. China

Service Hotline: 400-666-2522 Email: 3nh@3nh.com

Tel: 86-0755-26508999 Fax: 86-0755-27190609